Weight Scale with body composition analysis

User Manual

This scale allows you to assess your body composition. Knowing your body composition is helpful to combine proper nutrition and physical activity, which are essential for improving quality of life, good health and preventing disease.

-  **Read the instructions carefully before using this device.**

IMPORTANT INFORMATION ABOUT BODY COMPOSITION

It is important to know the body composition of people to understand their health status. Weight alone is not a representative indicator; thus, it is essential to know its composition. This means assessing two different aspects of body mass, as follows:

- Lean body mass: bones, muscles, intestines, water, organs, and blood
- Fat body mass: fat cells

Body fat

Total body fat mass is calculated by adding essential (primary) fat to stored (reserve) fat. Essential fat is needed for metabolism and a portion of stored fat, which in theory is in excess, is needed for physical activity. For men, a healthy body fat percentage is between 17-26% and for women between 22-31%. Body fat percentage increases with age.

Body water

The percentage of water in our bodies is an important indicator of a person’s general well-being, as water is the most important element in our bodies. Scientific research shows that 57% of water in the body is an optimal percentage, with a variation of +/- 10%. The percentage of water in the body is highest in children and decreases with age. Women, who typically have more fat tissue than men, have a body water percentage of approximately 55-58%, while men of approximately 60-62%. A percentage lower than 40% is considered too low, and above 70% too high. Changes in body weight within a few hours are considered normal and are due to water retention, as extra-cellular water is the only component that changes so quickly.

Muscle mass

Muscle mass is important in determining the body composition of a healthy person. A person with a high percentage of muscle mass will have less difficulty moving, but will also need more effort to do so. Physical exercise is essential for being fit, and the percentage of muscle mass is an effective indicator for controlling muscle mass. The percentage of muscle mass as part of the body weight is between 38-58% for men and 28-39% for women, depending on age and physical activity.

Bone mass

Bone mass increases quickly during childhood and reaches its maximum peak between the ages of 30 and 40. It then decreases slightly with ageing. A proper diet and regular physical exercise that includes also increasing muscle mass can help to reduce the bone degeneration processes. It is difficult to significantly affect bone mass. Factors such as age, gender, weight and height have little effect on bone mass. The average adult bone mass is 15% for men and 12% for women. There are no reliable guidelines or recommendations on bone mass.

Body mass index (BMI)

The body mass index is a weight-for-height index commonly used to classify underweight, overweight and obesity in adults. It is calculated by dividing the body weight (in kilograms) by the square of a person’s height (in metres). The World Health Organisation’s classification of the body mass index is as follows:

- Normal weight: 18.50-24.99 Kg/m²
- Underweight: less than 18.50 Kg/m²
- Overweight: over 24.99 Kg/m²
- Obese: over 29.99 Kg/m²

Basal metabolic rate (BMR)

The basal metabolic rate expressed in kilo-calories is the amount of energy required by the body,

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in fully resting conditions, to maintain its vital functions.

About 70% of the body’s energy consumption is due to vital organ processes. Approximately 20% of energy consumption comes from physical activity and a further 10% from the digestive process. It is possible to calculate this through an equation, using gender, age, height and weight. There are no reliable guidelines or recommendations on basal metabolic rate (BMR).

RECOMMENDATIONS FOR USE

- Place the scale on a flat, horizontal surface. Uneven or soft floors (e.g. carpets, linoleum) may cause inaccurate weight measurements.
- Step on the scale with your feet in a parallel position and distribute your weight equally.
- Do not move while the scale is measuring.
- It is important to weigh yourself with dry, bare feet, correctly positioned on the electrodes. To weigh yourself correctly, avoid touching your legs.
- The scale can be used as a standard scale by anyone. The following categories of people should not perform diagnostic weighing:
 - Pregnant women
 - People undergoing cardiovascular treatment
 - People with pace-makers or other electrical medical devices.
 - People running a fever
 - Children under 10 years of age
 - People subject to dialysis, with oedema of the legs or suffering from dimorphism or osteoporosis
- Always weigh yourself on the same scale, every day at the same time, preferably undressed and before breakfast. For best results, weigh yourself twice and if the two weights are different, take the average of the two readings. In the morning after getting up, wait at least 15 minutes before taking your weight measurement, especially measurements regarding body fat and body water in order to allow the water to distribute evenly throughout the body.
- The scale surface may be slippery if wet.
- The accuracy of the results may be compromised in the following cases:
 - presence of a high level of alcohol/drugs/caffeine in the body
 - after intensive fitness activity
 - when ill or pregnant
 - after a heavy meal or when dehydrated
 - if the personal data (height, age, gender) were not entered correctly
 - with wet or dirty feet

 **WARNING:** The body composition measured using this scale must only be considered as general and approximate information. For more details, contact your physician or dietitian.

 **WARNING:** The body composition measurement function of this scale can only be used by people between 10-100 years of age and 100-200 cm tall.

USING THE DEVICE FOR THE FIRST TIME

How to install the batteries

This scale operates with two 1.5V AAA batteries. When the batteries are low, the display shows “LO” (5). Replace the batteries in the appropriate compartment (4) and install them paying attention to the polarity.

Selecting unit of measure Kg-St-Lb

Select the desired unit of measure Kg, St or Lb when the scale is turned on and “0.0” is displayed. Press the ►▼ button (3) to change the unit of measure (otherwise it will not be possible to change it).

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Using the scale for weight measurement

- To obtain accurate measurements, place the scale on a firm surface.
- Step on the scale and wait without moving; after a few seconds, the display will show your weight (7). During the measurement, do not lean on anything to avoid inaccurate results.
- Step off the scale. The display will continue to show the weight for a few seconds, then it will turn off automatically.
- When the “ERR” symbol (18) appears, an error occurred during the measurement. Start the whole procedure again from the beginning.

ENTERING PERSONAL DATA

The scale can store the data of 10 users.

- Turn on the scale by pressing the SET button (1).
- «P1» (8) will flash on the display specifying user “1”: to change the memory location (from 1-10) press the ►▼ (3) button or ◄▲ (2) button. After selecting the chosen memory number, press SET (1) to confirm.
- Press the ►▼ (3) button or ◄▲ (2) button to confirm the gender (9). Press SET (1) to confirm.
- Use the ►▼ (3) button or ◄▲ (2) button to enter the age (from 10 to 100 years) (10). Press SET (1) to confirm.
- Enter the height (from 100-200 cm) using the ►▼ (3) button or ◄▲ (2) button to increase or decrease the values. Press SET (1) to confirm. The display will now show the stored data in sequence for verification. The scale is now ready for use.

To change the data entered, recall the related memory number and repeat the whole procedure. To enter data for another user, select an empty memory number and repeat the above procedure.

 **WARNING:** The scale automatically turns off if no button is pressed after 10 seconds.

CALCULATION OF BODY COMPOSITION

- Turn on the scale by pressing the SET button.
- «P1» (8) will flash on the display specifying user “1”. If you wish to change your memory storage number (from 1 to 10) press the ►▼ (3) button or ◄▲ (2) button. Select the chosen memory storage number, press SET (1) to confirm. The data entered and the values of the last measurements will be shown on the display. “0.0” (6) will then appear. To start the measurement immediately without displaying the data of the previous measurement, press SET (1) again and “0.0” will appear right away.
- Step on the scale placing your bare feet on the two electrodes. The display will show the following data in sequence:
 - the weight (7) flashes 3 times
 - the percentage of fat (12)
 - the percentage of water (13)
 - the muscle mass in Kg (14)
 - the percentage of bone mass (15)
 - the BMI (body mass index) in Kg/m² (16)
 - the BMR (basal metabolic rate) in Kcal (17)
- The scale will turn off automatically after about 5 seconds.

 **WARNING:** If the “ERR%” (19) symbol is displayed, start the whole procedure again from the beginning. When the “ERL” (20) symbol is displayed the body fat percentage is lower than the minimum measurable threshold. When the “ERL” (21) symbol is displayed the body fat percentage is higher than the maximum measurable threshold.

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